




0050



0050, Mini Desiccator Cabinet

Specifications:

- Injection molded PS cabinet is sturdy and unbreakable.
- Sides and door are completely transparent for easy viewing.
- The door AS made has a foam rubber gasket and roller latches for an airtight seal.
- When opened, the door slides into an out-of-the-way compartment beneath the cabinet.
- Smoky transparent cabinet is excellent for protecting light-sensitive materials.
- Suitable for the storage of highly sensitive components.
- Stackable to save bench top space.
- Materials: PS (Polystyrene)
- Compartment inside DIM.: W195 x D154 x H113mm
- (Accessories): 1 jar of moisture-absorbing silica gel in each compartment.

Model	Specifications	Outside dimension (mm)
0050	1 Level	W225 x D208 x H170



0040

0030/0040, Dry Keepers

Auto type adopts a high performance dehumidification unit ED-1 which maintains the Inside humidity at 30% to 40%RH. (With silica gel regeneration function. Troublesome replacement is unnecessary.) Non-auto type has a pan for silica gel and is convenient. Also use in applications other than dehumidification.

- Since there may be




0030

cases when the error produced by sudden humidity changes and shock becomes large, use the hygrometer supplied with the auto type only as a guide.

- The inside humidity is slightly different depending on the room temperature and usage conditions. Frame black. body smoky transparent.

Functions and strength are improved markedly by combining molded parts. Shelf plates can be adjusted in 15 steps at 28mm intervals to match the size of the stored items.

0030 Dimensions



Model	Type
0030	Latch door type Shelf plate 3 plates
0040	Latch door type Shelf plate 2 plates

0040 Dimension

