



SK310-230

SK310-230/SK410-230, Heat Sealer, sealing width – 300/400mm

- It fits nearly everywhere for world-wide use in laboratory.
- Polyethylene, polypropylene films can be easily sealed.
- Easy to handle with convenient and reasonable structure.
- Powerful transformer and high quality circuit components guarantee years operation.

Model	SK310-230	SK410-230	SK510-230
Available resin	polyethylene, polypropylene		
Sealing width	2mm		
Sealing length	300mm	400mm	500mm
Weight	4.3kg	5.4kg	6.5kg
Dimensions (mm)	W440xD85xH175	W530xD85xH175	W630xD85xH175
Wattage	390W	750W	900W
Power	AC 110V/220V, 50/60Hz		